

AIRTRAQ SP

Video laryngoscopy for remote locations and emergencies

Open and Intubate.
Fully Disposable



DESCRIPTION SIZE AND CODE	ET TUBE SIZES	MOUTH OPENING	COLOUR
REGULAR Size 3 A-011/ATQ-011	7.0-8.5	16 mm	 Blue
SMALL Size 2 A-021/ATQ-021	6.0-7.5	15 mm	 Green
PAEDIATRIC Size 1 A-031/ATQ-031	4.0-5.5	12 mm	 Purple
INFANT Size 0 A-041/ATQ-041	2.5-3.5	11 mm	 Grey
NASO TRACHEAL A-061/ATQ-061	Not applicable	15 mm	 Orange
DOUBLE LUMEN Endobronchial tubes A-071/ATQ-071	28-41 FR Left and right standadr and hooked	18 mm	 Yellow

AIRTRAQ AVANT Video laryngoscopy for routine intubations

Minimum cost per intubation.
Disposable Blades + Reusable Optics



Reusable optics

- **Airtraq Avant Optic (A-501)**
- Contains the optical, anti-fog and electronic systems.
- **When the battery is fully charged the Optics can be used for approximately 120 minutes.**
- The Optics include **light indicators** for **battery charge** and **service life**.



Low Cost disposable Blade in two different sizes



- **Disposable Blade and Eyecup**
 - Fully protects Optics from patient contact.
 - Quick and easy assembly with Optics.
 - Include lenses at the distal end of blade and also eyecup.
 - Sizes:
 - **Airtraq Avant Regular Blade : Size 3, Blue (A-511)**
For use with ETT 7.0 – 8.5. Minimum mouth opening: 17 mm
 - **Airtraq Avant Small Blade : Size 2, Green (A-521)**
For use with ETT 6.0 – 7.5. Minimum mouth opening: 17 mm

Optional accessories

Phone Adapter

- Universal Phone adapter (A-307)
- Install Free “Airtraq Mobile” App into iPhone® & Android Phor



Docking Station

- **Airtraq Avant Docking Station (A-590)**
- Recharges the battery.
- Displays service life of the optics.
- Can also be used to charge the Airtraq Wi-Fi Cameras

